

## IGH

Test Name	IGH
Also Known As	IGH
Panel?	No
Diagnosis Under Consideration	B-Cell Lymphoma

### Ordering Details

Test Code	5032-8
Assay Schedule	Mon - Sat, daily
Turnaround Time	7-10 days Turnaround
Preferred Specimen and Volume	2 mL Heparinized or EDTA Marrow, 5 mL Heparinized Blood or Fresh Tissue in RPMI.
Alternate Specimen and Volume	Non- decalcified Paraffin Embedded Tissue ( block)

### Specimen Handling Information

Specimen Storage	Refrigerate fresh specimen until shipping, ship within 24 hrs.
Specimen Transport	Ship specimen with ice pack.
Causes For Specimen Rejection	No patient ID on Specimen, Patient ID on specimen and Requisition doesn't match, incorrect fixative (tube/solution), Specimen without GenPath/BioReference requisition. Frozen specimen.

### Billing Information

CPT Codes	88368 (1)
-----------	-----------