

X/Y

| | |
|-------------------------------------|---|
| Test Name | X/Y |
| Also Known As | X/Y |
| Panel? | Yes |
| Clinical Utility and/or Limitations | Used to monitor Bone Marrow Transplants |
| Diagnosis Under Consideration | BM Transplant Monitoring |

Ordering Details

| | |
|-------------------------------|---|
| Test Code | 5066-6 |
| Assay Schedule | Mon - Sat, daily |
| Turnaround Time | 24-72hr Turnaround |
| Preferred Specimen and Volume | 2 mL Heparinized or EDTA Marrow, 5 mL Heparinized or EDTA Blood or Fresh Tissue in RPMI, or sorted cells in sterile tube. |
| Alternate Specimen and Volume | Non- decalcified Paraffin Embedded Tissue (block) |

Specimen Handling Information

| | |
|-------------------------------|--|
| Specimen Storage | Refrigerate fresh specimen until shipping, ship within 24 hrs. |
| Specimen Transport | Ship specimen with ice pack. |
| Causes For Specimen Rejection | No patient ID on Specimen, Patient ID on specimen and Requisition doesn't match, incorrect fixative (tube/solution), Specimen without GenPath/BioReference requisition. Frozen specimen. |

Billing Information

| | |
|-----------|-----------|
| CPT Codes | 88368 (1) |
|-----------|-----------|