

Anti-Phosphatidyl Ethanolamine Antibodies (IgG, IgM)

Test Name	Anti-Phosphatidyl Ethanolamine Antibodies (IgG, IgM)
Panel?	Yes
Tests in Panel	IgG, IgM

Ordering Details

Assay Schedule	Once a week
Turnaround Time	5-7 Days
Preferred Specimen and Volume	1 ML of Aliquoted Frozen Serum (from the Red Top Tubes)

Specimen Handling Information

Specimen Storage	Keep Frozen
Specimen Transport	Ship Frozen sample packed on Dry Ice by a reliable courier.
Causes For Specimen Rejection	Specimen will be rejected if grossly hemolyzed.

Billing Information

CPT Codes	83520 (2)
-----------	-----------