

## Dilute Russell Viper Venom Test- Complete

Test Name	Dilute Russell Viper Venom Test- Complete
Also Known As	dRVVT - Complete
Panel?	No

### Ordering Details

Assay Schedule	Monday- Friday, Daily
Turnaround Time	1-3 days
Preferred Specimen and Volume	1 ML of Aliquoted Frozen Plasma (from the Blue 3.2% Citrate Tubes). Note: Blue Citrate Tubes must be filled to proper level before spinning / aliquot.

### Specimen Handling Information

Specimen Storage	Keep Frozen
Specimen Transport	Ship Frozen sample packed on Dry Ice by a reliable courier.
Causes For Specimen Rejection	Specimen will be rejected if grossly hemolyzed, clotted, or thawed.

### Billing Information

CPT Codes	85613 (2)
-----------	-----------

### Other Information

Notes	Screen & confirmatory reagent.
-------	--------------------------------