

D-Dimer, Quantitative

Test Name D-Dimer, Quantitative

Ordering Details

Assay Schedule Monday- Friday, Daily

Turnaround Time 1-3 days

Preferred Specimen and Volume 1 ML of Aliquoted Frozen Plasma (from the Blue 3.2% Citrate Tubes). Note: Blue Citrate Tubes must be filled to proper level before spinning / aliquot.

Specimen Handling Information

Specimen Storage Keep Frozen

Specimen Transport Ship Frozen sample packed on Dry Ice by a reliable courier.

Causes For Specimen Rejection Specimen will be rejected if grossly hemolyzed, clotted, or thawed.

Billing Information

CPT Codes 85379