

## Factor IX (9) Clotting Activity

Test Name Factor IX (9) Clotting Activity  
Also Known As Christmas factor

### Ordering Details

Assay Schedule Monday- Friday, Daily  
Turnaround Time 2-5 days  
Preferred Specimen and Volume 1 ML of Aliquoted Frozen Plasma (from the Blue 3.2% Citrate Tubes). Note: Blue Citrate Tubes must be filled to proper level before spinning / aliquot.

### Specimen Handling Information

Specimen Storage Keep Frozen  
Specimen Transport Ship Frozen sample packed on Dry Ice by a reliable courier.  
Causes For Specimen Rejection Specimen will be rejected if grossly hemolyzed, clotted, or thawed.

### Billing Information

CPT Codes 85250