

Abnormal PT & aPTT Assessment

Test Name	Abnormal PT & aPTT Assessment
Panel?	Yes
Tests in Panel	5798- PT, 5835- aPTT, 5830- aPTT mixing study, 5840-PT mixing study, 5725- Factor V (5) Clotting Activity,5735- Factor X (10) Clotting activity,5743- D-Dimer “ Quantitative, 5724- Factor II (2) Clotting Activity, 5728- Factor VII (7) Clotting Activity,5745- Fibrinogen Clottable Activity,5807- Thrombin Time (TT), 5729- Factor VIII (8) Clotting activity, 5769- Factor VIII (9) Clotting activity, 5737- Factor XI (11) Clotting activity,5771- Partial Thromboplastin Time (aPTT) Lupus Sensitive, 5730-dRVTT- Complete, 5749- Hexagonal Phase Assay

Ordering Details

Assay Schedule	See individual tests
Turnaround Time	5-7 days
Preferred Specimen and Volume	Draw 5 blue 3.2% citrate tubes to submit 6 aliquot Frozen Plasma (FP) tubes.

Specimen Handling Information

Specimen Storage	Frozen Plasma aliquot tubes must be frozen and shipped on dry ice.
Causes For Specimen Rejection	Frozen Specimens will be rejected if grossly hemolyzed, clotted, or thawed.

Billing Information

CPT Codes	85732 (9), 85611 (8), 85230 (1), 85610 (1), 85730 (1), 85220 (1), 85260 (1), 85384 (1), 85379(1), 85210 (1), 85670 (1), 85240(1), 85250(1), 85270(1), 85613 (2), 85597 (2)
-----------	--