

Antiphospholipid Antibody and Lupus Inhibitor Assessment- Expanded

Test Name	Antiphospholipid Antibody and Lupus Inhibitor Assessment- Expanded
Panel?	Yes
Tests in Panel	5780- CBC - Diff & Platelets,5798- PT, 5835- aPTT,5710- Anti-Beta2 Glycoprotein I Abs (IgG, IgM, IgA), 5720- Anti-Cardiolipin Abs (IgG, IgM, IgA), 5790- Anti-Phosphatidyl Ethanolamine Abs (IgG, IgM)5820- Anti-Prothrombin Abs (IgG,IgM), 5706- Anti-Annexin V (5) Abs (IgG, IgM), 5790- Anti-Phosphatidylserine Abs (IgG, IgM)Mixing Studies - 5840-PT or 5830- aPTT (If Protome prolonged,add 5724- Factor II (2) clotting activity)5730- dRVVT - Complete, 5749- Hexagonal Phase Assay,5771- aPTT Lupus Sensitive, 5784-Platelet Neutralizing Procedure

Ordering Details

Turnaround Time	7-10 days
Preferred Specimen and Volume	Draw 1 lav, 2 red top tubes and 3 blue 3.2% citrate tubes to submit 1 Lav (room temp) 4 Frozen Serum (FS) aliquot tubes and 4 Frozen Plasma (FP) aliquot tubes.

Specimen Handling Information

Specimen Storage	Lavender Tubes should be stored and shipped separately at Room Temperature. Frozen Plasma and Frozen Serum aliquot tubes must be frozen and shipped on dry ice.
Causes For Specimen Rejection	Frozen Specimens will be rejected if grossly hemolyzed, clotted, or thawed.

Billing Information

CPT Codes	85025, 85610, 85732 (9), 86147 (3), 86146 (3), 85611 (8), 85730, 85597 (3), 86148 (2), 85613 (2),0030T(2),83520(2),83516(2)
-----------	---