

DIC Assessment

Test Name	DIC Assessment
Panel?	Yes
Tests in Panel	5780- CBC - Diff & Platelets,5798- PT, 5835- aPTT,5824- Fibrin Monomer,5745- Fibrinogen Clottable Activity,5743- D-Dimer, Quantitative,5797- Prothrombin Fragment 1.2 (F1.2),5753- Smear Review (for morphology),5808- Thrombin/Anti-Thrombin Complex (TAT),5807- Thrombin Time (TT),

Ordering Details

Turnaround Time	15 Days
Preferred Specimen and Volume	Draw 1 lavender and 3 blue 3.2% citrate tubes to submit 1 Lav (room temp) and 4 Frozen Plasma (FP) aliquot tubes.

Specimen Handling Information

Specimen Storage	Lavender Tubes should be stored and shipped separately at Room Temperature. Frozen Plasma aliquot tubes must be frozen and shipped on dry ice.
Causes For Specimen Rejection	Frozen Specimens will be rejected if grossly hemolyzed, clotted, or thawed.

Billing Information

CPT Codes	85025, 85610, 85730, 85384, 85379, 85670, 85362, 85060, 83520(2)
-----------	--