

ER/PR/Ki-67, Plus DNA (Technical Only)

Test Name	ER/PR/Ki-67, Plus DNA (Technical Only)
Also Known As	ER/PR/Ki-67, Plus DNA (Technical Only)
Panel?	Yes
Clinical Utility and/or Limitations	Breast Cancer Treatment Stratification & prognosis
Diagnosis Under Consideration	Breast Cancer

Ordering Details

Test Code	5129-2
Assay Schedule	Monday-Saturday
Turnaround Time	24 Hrs.
Preferred Specimen and Volume	Paraffin Embedded Tissue
Alternate Specimen and Volume	4-7 micron thick unstained positively charged slides

Specimen Handling Information

Specimen Storage	Room Temperature
Specimen Transport	Ship specimen ambient except during the months of May to September ship with ice pack to prevent damage from the heat.
Causes For Specimen Rejection	No Tumor in Block. No Patient Id and/or Surgical number on Specimen, Patient ID and/or Surgical number on specimen and Requisition doesn't match, incorrect fixative (tube/solution), Specimen without GenPath/BioReference requisition.

Billing Information

CPT Codes	88360 TC x 3, 88182 x 1
-----------	-------------------------