

ER/PR/Ki-67, Plus DNA (Technical only - available only to Pathologists)

Test Name	ER/PR/Ki-67, Plus DNA (Technical only - available only to Pathologists)
Also Known As	ER/PR/Ki-67/HER2 - TC , Plus DNA
Panel?	Yes
Clinical Utility and/or Limitations	Breast Cancer Treatment Stratification and Prognosis
Diagnosis Under Consideration	Breast Cancer

Ordering Details

Test Code	5131-8
Assay Schedule	Monday-Saturday
Turnaround Time	24 Hrs.
Preferred Specimen and Volume	Paraffin Embedded Tissue
Alternate Specimen and Volume	4- 7 micron thick unstained sections on positively charged slides and 60 micron shavings, 3 pieces

Specimen Handling Information

Specimen Storage	Room Temperature
Specimen Transport	Ship specimen ambient except during the months of May to September ship with ice pack to prevent damage from the heat.
Causes For Specimen Rejection	No Tumor in Block. No Patient Id and/or Surgical number on Specimen, Patient ID and/or Surgical number on specimen and Requisition doesn't match, incorrect fixative (tube/solution), Specimen without GenPath/BioReference requisition.

Billing Information

CPT Codes	88182 x 1, 88360-TC x4
-----------	------------------------