

DNA Ploidy for Molar Pregnancy by FLOW Cytometry

Test Name DNA Ploidy for Molar Pregnancy by FLOW Cytometry
Diagnosis Under Consideration Assess for Partial Hydatidiform Mole

Ordering Details

Test Code 5555-8
Assay Schedule Monday-Saturday
Turnaround Time 24 Hrs.
Preferred Specimen and Volume Paraffin Embedded Tissue
Alternate Specimen and Volume Fresh Tissue or 60 micron shavings (4 pieces)

Specimen Handling Information

Specimen Storage Room Temperature
Specimen Transport Ship specimen ambient except during the months of May to September ship with ice pack to prevent damage from the heat.
Causes For Specimen Rejection No Tumor in Block. No Patient Id and/or Surgical number on Specimen, Patient ID and/or Surgical number on specimen and Requisition doesn't match, incorrect fixative (tube/solution), Specimen without GenPath/BioReference requisition.

Billing Information

CPT Codes 88182 x 1