

Arterial and Atherothrombotic Risk- Basic

Test Name	Arterial and Atherothrombotic Risk- Basic
Panel?	Yes
Tests in Panel	Lipid profile,5755- Homocysteine,5719- hs C-reactive protein,5758- Lipoprotein a- Lp(a),5745- Fibrinogen Clottable Activity,5798- PT,5807- Thrombin Time (TT),5835- aPTT,5729- Factor VIII (8) clotting activity,
Clinical Utility and/or Limitations	Assessment of patients who have had an arterial thrombosis or family history of early onset atherothrombotic disease
Diagnosis Under Consideration	Panel for patients with stroke, coronary artery disease or peripheral artery disease.

Ordering Details

Test Code	5874-3
Turnaround Time	5-7 days
Preferred Specimen and Volume	Draw 2 red top tubes and 2 blue 3.2% citrate tubes to submit 4 Frozen Serum (FS) aliquot tubes and 6 Frozen Plasma (FP) aliquot tubes.

Specimen Handling Information

Specimen Storage	Frozen Plasma and Frozen Serum aliquot tubes must be frozen and shipped on dry ice.
Causes For Specimen Rejection	Frozen Specimens will be rejected if grossly hemolyzed, clotted, or thawed.

Billing Information

CPT Codes	83090, 83695, 85384, 86141, 85610, 85730, 85240, 80061, 85670
-----------	---------------------------------------------------------------